Precise and reliable –
we know what is
important to you.

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Poulten & Graf Volumetric Glassware - tested precision-made glassware

Poulten & Graf offers a wide range of Volumetric Glassware under the trade marks FORTUNA® and .

The whole range of Volumetric Glassware made by Poulten & Graf is produced according to British BS, European EN or international ISO standards.

There is also a specialised range made according to US Pharmacopoeia (USP) standards. In addition  cylinders and flasks are made to ASTM standards.

Poulten & Graf Volumetric Glassware is produced according to the highest quality standards. Control of raw materials, production and quality, work hand in hand and guarantee the high quality of the final product.

Poulten & Graf Volumetric Glassware, conformity certified, class A, is produced with a batch number and can be delivered with a batch certificate (except for Original FORTUNA® graduated and bulb pipettes with suction piston are delivered without certificate).

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Graduation

We distinguish between two different graduation types:

- for Class A or AS products: complete ring marks
- for Class B products: short graduations

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Calibration

Generally we differentiate between volumetric glassware which is calibrated to contain (IN) or calibrated to deliver (EX). Flasks, measuring and mixing cylinders are calibrated to contain (IN) so the liquid poured out corresponds with the volume indicated. Pipettes and burettes, however, are calibrated to deliver (EX). Provided the required waiting time is observed, the delivered liquid corresponds with the volume indication.

FORTUNA® and  glassware is individually calibrated (no mass production technique and no approximations).

Special Note: When users wish to check the accuracy of volumetric glassware, care should be taken not only to avoid parallax error, but also to avoid a false view of the meniscus. Use of a black ring or a clip approx. 1 mm below the mark will give a black reflection on the bottom of the meniscus and will dramatically improve the accuracy of the reading (see ISO 4787).

We can supply suitable black clips on request.

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Colour marking

The colour markings on FORTUNA® and **volac** volumetric glassware are applied by silk screen printing. We only use diffusion or enamel colours. The diffusion colour - using amber stain - is completely absorbed into the glass surface during the firing process. It cannot be removed except by processes which remove the glass surface itself. It is therefore resistant to cleansing and even to very aggressive chemicals.

All other colours - e.g. blue, white, yellow, green - are enamel or melting colours, which are melted onto the glass surface at high temperature.

volac glassware is individually printed - both Class B and Class A. Particularly for graduated glass (graduated pipettes, burettes, cylinders) individual printing is the only way to ensure that the top and bottom of the scale matches individual calibration exactly.

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Classes

In accordance with the ISO standard: FORTUNA®

Class	Certificate	Individually Numbered	Products
A	batch certificate	on request	cylinders, flasks
AS	batch certificate	on request	burettes, graduated and bulb pipettes

In accordance with the ASTM, ISO and USP standards: 

Class	Certificate	Individually Numbered	Products
B	no	no	burettes, graduated and bulb pipettes, cylinders
A	no	with or without individual certificate	burettes, graduated and bulb pipettes, cylinders, flasks
ASTM	single certificate	yes	cylinders, flasks
USP	single certificate	yes	burettes, graduated and bulb pipettes, cylinders, flasks

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Why use volumetric glassware according to ASTM- or USP- standard?

In the United States and in Arabic Countries, as well as in many Asian Countries there is strong demand for volumetric glassware according to ASTM- or USP-standards.


Essentially, the tolerances of measuring devices made to these standards differ to a great extent from the European norm, ISO, and are much tighter.

If a company produces goods that have to comply with the American FDA import controls or come within their responsibility, they have to meet certain regulations (GMP).

GMP stands for Good Manufacturing Practice. This is a set of regulations, which ensures good quality, good production processes and appropriate production environments in the production of pharmaceutical products and chemical agents, as well as cosmetics, food and animal feed.

Amongst other things the focus lies in a complete and closed traceability of the products. In addition, the traceability and control of the manufacturing process (tracking), plays an important role.

To ensure the tracking of the volumetric glassware which Poulten & Graf produces, individual serial numbers and quality certificates are provided. Manufacturers, who do not follow the GMP and cannot ensure complete quality control, face warnings and penalties. As a last resort, producers, who aim to sell their products mainly in territory controlled by the FDA, can lose their permission to sell into the US market. To ensure maximum reliability, precision and traceability, choose **volac** volumetric glassware, produced according to ASTM- or USP- standard, with individual serial numbers and certificates.

volac Volumetric Glassware Specifications

Volumetric Flasks Class A	ASTM and USP is tighter than ISO
Bulb Pipettes Class A	some USP are the same as ISO but some are tighter
Graduated Pipettes Class A	again some USP specifications are tighter than ISO but some are the same
Measuring Cylinders Class A	ASTM and USP specifications are unique
Burettes Class A	USP and ISO are identical

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