Certificate of Registration

nqa

This is to certify that the Quality Management System of

Proton & Seaf Limited
Plot 10/10A, T. Alakshy Day Industrial Estate, Alkhuw, Barking, Essex, IG11 0AG,
available to

Manufacture and sale of laboratory glassware and equipment

has been assessed and registered by NQA against the provisions of

BS EN ISO 9001 : 2008

The registration is subject to the company maintaining a quality management system,
to the above standard, which will be monitored by NQA.

UKAS

Certification Directory

UKAS

Certificate No.
Date
Registered
NQA 1045
Date Issue

2012
7 September 2012
24 June 2013
27 June 2015
2018

This certificate is issued on the basis of the information provided by the client and is not intended to be used for any other purpose. It is not a guarantee of performance and does not constitute an endorsement of the client or its products. The certificate is issued on the basis of the information provided by the client and is not intended to be used for any other purpose. It is not a guarantee of performance and does not constitute an endorsement of the client or its products.

Poulsen & Graf GmbH

What is understood by quality?

In general usage the term 'Quality' has a number of meanings. Glass products of first class quality might be called this because the finish is good, the raw material used is good, the packaging is good, the appearance is good or the accuracy is good. Quality is defined as having all those characteristics and features of a product, which refers to its suitability for the fulfilment of given requirements, that is 'fit for purpose'. So Quality is related to the need and the specification set.

The Quality Management of Poulten & Graf

Product quality is directly dependent on the quality of process. Quality assurance must begin in the very first phase of the product creation process and cover all following phases. This global view of quality assurance, involving all aspects of production, design, testing, packing and conformity with both internal standards and externally set specifications, is the basis of our quality management. This is reflected in the approvals of DIN EN ISO 9001/2008, which we have achieved and which we maintain.

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06. Technical Support · Technische Details

Certification or attestation of conformance

Since January 01, 2015, new German Metrology and Verification Regulations have come into force. They replace the previous German Weights and Measures Regulations. For Class A and Class AS Volumetric Glassware which are produced in Germany, a new conformity sign is being used. The H-sign is now replaced by the **DE-M**-sign.

As before, the manufacturer is allowed to certify the conformity on basis of the PTB-regulations (National Metrology Institute). Government institutions will ensure that the regulations are being followed. Therefore, nothing has changed in the production and testing methods. If, required, a manufacturer's declaration of conformity can be enclosed with class A and AS volumetric products.

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Works Certification of Accuracy of Volumetric Glassware

All  Class A glassware is always individually tested. Test records are kept in house with the manufacturing records. However individual Works Certificates of Accuracy can be provided for each piece if required; each piece of glass has its own individual number engraved on it and the Certificate supplied will have the corresponding number, the date and the relevant test results.

Such Works Certificates of Accuracy (look for the designation WAC in the product codes) should not be confused with Certificates of Conformity which refer to batch testing only.

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The system for lot numbers on FORTUNA® volumetric glassware

All FORTUNA® conformity approved volumetric glassware is delivered with the FORTUNA® Batch Certificate of Performance. A regular control of all testing equipment is required by DIN EN ISO 9001. The requirement for traceability is met by the batch certificate and with the easily identifiable lot numbers (e.g. LOT B008) on FORTUNA® volumetric glassware, which is produced in Wertheim.

Beside the conformity sign all FORTUNA® volumetric glassware is marked with an individual lot number that allows an exact classification back to the batch certificate.

If your certificate is lost, you may ask for a copy by informing us of the product details and the batch number or request it online.

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R0009

To define the different classes (Class A, Class AS, Class B, etc.) for glass or 'performance' for pipet-tors and dispensers the following terms are used: 'accuracy' and 'coefficient of variation' (precision). Sometimes we also speak about 'tolerances' which are speci-fied in the ISO-standards. This means the permitted deviation of the measured values from the nominal volume. By accuracy we mean corresponding to the aver-age value of a series of measure-ments adjusted for temperature and pressure with the desired value. Precision is a measure of the variation from the nominal value achieved through a series of measurements at one setting. The coefficient of variation is the % value of the precision.

Measurement of volumes is part of the routine work in a laboratory. The most com-mon volumetric glassware such as bulb- and graduated pipettes, burettes, volumetric flasks and graduated cylinders form part of the basic equipment in a laboratory. Inscriptions and graduations on volumetric glassware should be easy to read and acid- and alkali-resistant. The complete range of FORTUNA and  volumetric glassware fulfils these requirements.

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
Calibration of Volumetric Glassware

- **EX** calibrated to deliver. The delivered liquid corresponds with the volume indication. The liquid that remains on the wall or in the tip of the glassware is taken into consideration.
- **EX + s** calibrated to deliver after a waiting time (for pipettes 5 seconds, for burettes 30 seconds). Keeping to the waiting time is important to avoid mistakes in measurement.
- **IN** calibrated to contain. The contained liquid corresponds with the volume indication. Example: flasks and cylinders.
- **20°C** The calibration has been done at 20°C. Volumetric glassware changes volume about 0,1 % each degree Celsius.

Precision Classes

- **A** tolerance according to ISO standard. Class A is the most precise class. Glassware of class A is suitable for official calibration, this means it can be tested by either the German or the UK Weights and Measures Offices. This is useful for control of measuring instruments according to DIN ISO 9001.
- **AS** as class A, but for pipettes and burettes with fast delivery. The waiting time is much shorter than for class A.
- **B** the tolerances of class B are about twice those of class A and AS. Volumetric glassware of class B is suitable for routine work in the laboratory.

Meniscus

Value is read at the lowest point of the liquid surface, which means when the lowest point of the meniscus touches the upper line of the graduation mark (see picture 1).

For burettes with Schellbach stripes the reading has to be made at the level where the two wedge-shaped points are touching (see picture 2). When reading the scales it is important that one's eyes are at height of the fluid level. Otherwise it results in reading errors (parallax errors).

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
A Guide to calibrate Micro Pipettors and Bottle Top Dispensers (compiled from BS ISO 8655: Part 1-6)

Test Lab Requirement

- Lab should ideally be air conditioned.
- Temperature and air pressure should be stable.
- Air movement should be minimised.
- Weighing/balance surface should be free of vibration.

Test Lab Equipment

- Air pressure barometer
- Certified Class A thermometer
- Accurate balance (5 decimal places)
properly calibrated to National Standards
- Distilled water at stable room temperature

Pipetting Practice for Accurate Results

Use high quality tips preferably of a brand to match the pipettor. Ensure tips do not leak (hold pipettor vertically, tip filled with liquid, without dripping for at least 20 secs). Rinse and pre-wet each tip with distilled water at least 3 times before use. Fill and dispense slowly and smoothly. Fill with pipettor tip not more than 3-4 mm below liquid surface and pipettor held vertical (fig. 1) Dispense with pipettor tip against wall of container and pipettor held at 45° to the vertical and draw tip up the container wall one or two seconds after second blow out (fig. 2).

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

06. Technical Support · Technische Details

Pipetting - Balance Readings

Tare a container on the balance. Using a wetted tip, pipette/dispense into the tared container on the balance. Read weight of dispensed water. Read temperature of distilled water. Read the air pressure. For setting calibration test 4 times at each volume, for checking/certifying calibration test minimum 10 times and for testing coefficient of variation test 30 times.

Calculation of Accuracy

- Calculate mean weight $\overline{W}$ from the number n of individual weighings W_i summed to give $\sum W_i$ using the equation:

$$\overline{W} = \frac{\sum W_i}{n}$$

- Calculate mean volume V_t of liquid samples at test temperature t :

$$\overline{V}_t = \overline{W} \cdot Z$$

where Z is the conversion factor for the density of water at given temperatures and air pressure - see table (page 294)
N.B. If an electronic balance is used, the effects of evaporation are minimal and can be ignored; if a mechanical balance is used mean evaporation loss $\overline{e}$ needs to be estimated and added to $\overline{W}$.

- Calculate inaccuracy in percent of the mean $\overline{A}_t$ at test temperature t using V_t mean volume and V_o nominal volume by equation:

$$A_t = \overline{V}_t - V_o \cdot \frac{100}{V_o}$$

Index of Symbols

$\overline{W}$	mean weight	$\overline{e}$	mean evaporation loss
W_i	results of individual weighings	Z	Conversion factor (see table)
n	number of weighings	$\overline{A}_t$	Accuracy (inacc.) as percent of mean
t	water temp. of test	S	Standard deviation of weighings (or volumes)
$\overline{V}_t$	mean volume at test temperature	CV	Coefficient of variation
V_o	nominal volume		

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06. Technical Support · Technische Details

Calculation of Precision (Repeatability) - Coefficient of Variation

- Calculate standard deviation **S** by the equation:

$$S = \sqrt{\frac{\sum (W_i - \overline{W})^2}{n-1}}$$

(n should be not less than 10
and should ideally be 30)

- Calculate coefficient of variation **CV**:

$$CV = \frac{S}{\overline{W}} \cdot 100$$

If the balance is attached to computing equipment $\overline{W}_i$ and **W** may be replaced by equivalent volumes.

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Calibration of Volumetric Glassware (see BS 6696)

Test Lab Requirement

- Full temperature control (air conditioning)
- Atmospheric Pressure Barometer
- Certified Thermometer
- Accurate calibrated weighing balance

Cleaning and Drying

Wash, clean and dry glassware thoroughly. Ensure that a satisfactory method of cleaning is employed, since any contamination will effect the meniscus shape/level.

During calibration it is critical that all internal surfaces are clean and dry to avoid moisture droplets or contaminants.

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
Timing Delivery

Delivery glassware should be timed before calibration, to ensure that total delivery time is within tolerances specified by the relevant BS and ISO standards.

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Weigh and Tare


With delivery glassware such as pipettes and burettes, weigh and tare a suitable sized weighing vessel.

Glassware should be maintained at a controlled room temperature (air conditioning).

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
Fill with purified water a few millimeters above the graduation mark. Pipettes and burettes should be filled from below, i.e. through the tip. If burettes are filled from the top for added convenience, take care to avoid wetting the walls above the upper graduation mark and check that no air is trapped in the stopcock. For pipettes and burettes control the flow of excess liquid from the tip to allow the meniscus to settle exactly on the graduation mark. For flasks and cylinders remove excess water to match meniscus and graduation mark exactly using a pipette and bulb.

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
When setting the meniscus ensure that the glassware is vertical. Ensure parallax error is avoided by viewing the meniscus from an eyelevel height against a white background.

In addition a black shade should be secured either around the article to be tested or placed behind it. The shade should be below the setting height (approx 1 mm lower). Ensure the meniscus edge is smooth, straight and level.

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06. Technical Support · Technische Details

For graduated/bulb pipettes and burettes, dispense water into a tared weighing vessel and weigh. For flasks and cylinders weigh the tared and filled glassware. Record the result from the balance. Use a certified Class A thermometer to measure the water temperature. Add the water temperature correction figure to the net weight obtained.

Measure the barometer reading and air temperature and add (or subtract if minus) the correction figure. Please ask us for the tables.

This produces the final gravimetric result which has already been converted to the volumetric figure in millilitres.

All of the procedure described must be followed to ensure that consistently accurate results are obtained.

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Environment and safety

Current discussions on environmental protection, dangerous materials etc. have made laboratory technicians more sensitive in the selection of laboratory instruments. Safety aspects remain the primary concern while handling dangerous or hazardous materials. Laboratory equipment must be extremely safe and reliable. All regulations must be observed to ensure protection of both humans and the environment.

We offer a repair service for a number of products:

- bottle top dispenser POLYFIX®, OPTIFIX®,
- electronic burette OPTIFIX® titrier,
- automatic dosing system OPTIMAT® 2,
- pipettors.

Please remember that when returning any instrument for repair it must be accompanied by a written certificate of decontamination.

When disposal of waste is involved we ask you to consider the necessary regulations. These apply also to disposal of single use products.

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06. Technical Support · Technische Details

Important safety instructions:

1. Laboratory instruments and equipment should be checked on purchase as being in good condition and working order before use.
2. Laboratory instruments should be checked for damage as part of a daily routine to avoid risk of injuries such as cuts, burns and infections.
3. Laboratory glassware which shows any sign of damage should be disposed of in a safe manner.
4. Sudden changes in temperature should be avoided while handling glassware.
5. Glass equipment should be inspected thoroughly before using in a vacuum or under pressure.
6. Sudden changes in pressure should be avoided. Vessels should be vented smoothly when used in a vacuum. Flat bottomed glassware should not be used under pressure or in a vacuum.

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06. Technical Support · Technische Details

Cleaning of Laboratory Glassware

Laboratory glassware should be washed before first use.

Cleaning by hand

The common cleaning method is to wipe and rub the glassware with a cloth or sponge soaked in a cleaning solution. Abrasive cleaners and sponges should not be used because they can damage the surface of the glass.

Machine cleaning

Cleaning in a machine is more gentle for glassware than soaking it. Glass comes into contact with the dishwashing liquid for relatively short periods when it is being sprayed onto the surface of the glass.

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Cleaning of Bottle Top Dispenser

To ensure proper operation and life of the OPTIFIX® dispensers, we recommend the following cleaning procedure to provide you with years of trouble free service.

- Prepare a standard cleaning solution and hot water.
- Remove the piston from the Cylinder via the volume setting knob. Remove discharge tube, intake tube and adapter.
- Put all parts into the cleaning solution for 2-4 hours or overnight to allow a good soaking.
- After soaking, use a clean glassware brush to remove any further dirt if necessary.
- Rinse all dispenser parts with distilled water and air dry.
- Reassemble the unit and put the OPTIFIX dispenser back into service.

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06. Technical Support · Technische Details

Disinfection

Laboratory equipment that has been in contact with infectious material should be disinfected before use, for the protection of laboratory workers. It can be washed by hand with a cleaner-disinfectant. Further it can be cleaned by using physical or thermal methods (for 10 minutes at 93°C as required by the BGA). If necessary the glassware can be steam sterilized.

Steam sterilization

Steam sterilization is 'the destruction or irreversible inactivation of all microorganisms of 120°C and 2 bar'. At a sterilization temperature of 121°C the glassware should be sterilized for 20 minutes minimum effective application time. Laboratory equipment should be always cleaned carefully before steam sterilization, otherwise any soiling will bake on during the steam sterilization. If there are any chemicals on the surface of the glassware during the sterilization they can damage the glass due to the high temperatures. Containers should be opened during the sterilization to avoid a build-up of pressure. The steam has to have unrestricted access to all contaminated points to ensure effective steam sterilization.

Raw materials we use

Soda Lime glass finds its use in the foodstuffs industry as bottle and preserve glass as well as in the laboratory industry. The chemical composition of this glass consists of 71-75 % sand (SiO_2), 12-16 % soda and 10-15 % lime. Soda lime glass is subject to acute stress under heat shock and has a relatively high coefficient of expansion. It is particularly resistant to alkali liquids.

Borosilicate glass differs from soda lime glass by using a higher proportion of sand (SiO_2). The composition of this glass contains 70-80% silicon dioxide (SiO_2), 7-13% boron trioxide (B_2O_3), 4-8% Na_2O and K_2O , as well as 2-7% alumina oxide (Al_2O_3). Glass with such composition is characterised by a high resistance against acidic attack, low thermal expansion and excellent temperature stability. This glass is particularly suitable for laboratories and chemical use where contact with aggressive materials is likely to be prolonged.

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06. Technical Support · Technische Details

The POLYFIX® Dispenser is made of Polypropylene and glass.

The OPTIFIX® dispenser is made of PTFE and glass, while the SMART pipettor is made of Polypropylene and Polycarbonate.

PP - Polypropylene

When polymerizing, propylene develops a product which contains methyl side groups of monomers of the paraffin chain. The volume requirement due to the regular methyl side chains is relatively large. For this reason polypropylene is the plastic with the lowest specific weight in compact form. It is similar to polyethylene, and may be used instead of polyethylene, where it particularly depends on good heat resistance, high impact strength, inherent stability during stress or particularly low specific weight. Polypropylene also has a good resistance to many inorganic reagents as well as the less aggressive organic solvents.

PTFE - Polytetrafluorethylene

PTFE has no thermoplastic characteristics although it consists of linear, chemical chain molecules which are not linked. It is not possible to handle or form PTFE in the usual thermoplastic way. Nevertheless PTFE is a thermoplast because it can be sintered. For the production of shaped parts from PTFE the polymer is cold pressed as powder under high pressure into the desired form and sintered at temperatures around 380 °C. PTFE is attacked only by elemental fluorine and chlorine tri fluoride at higher temperatures and pressures as well as by melting alkali metals. Otherwise it is resistant to all chemicals. A further advantage of polytetrafluorethylene: it has the lowest coefficient of friction of all solids. However above 350 °C PTFE burns in air to form noxious gases.

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