



Glass Disposable Pasteur Pipettes

Length

230 mm

unplugged,
heavy wall
soda lime glass

Art.-No.

D812/HW

Packing unit: 800



Length

150 mm

50-Dropper
unplugged
soda lime
glass

Art.-No.

D897

50-Dropper
pre-plugged
soda lime
glass

Art.-No.

D897/PP

Packing unit: 1000 (4 inner packs x 250)



Teat for Pasteur Pipettes

made of PVC, smooth for safety and durability

Art.-No.

D813



Packing unit: 100



Poulten & Graf Ltd

Peak Works, 1 Alfreds Way
Barking, Essex,
IG 11 0AS, UK
Tel.: +44 (0)20 8594 4256
Fax: +44 (0)20 8594 8419
www.poulten-graf.co.uk
VOLAC@poulten-graf.com

FORTUNA®

Poulten & Graf GmbH

Karl-Carstens-Str. 10,
97877 Wertheim, Germany
Tel.: +49 9342 9229-0
Fax: +49 9342 9229-80
www.poulten-graf.de
sales@poulten-graf.com
www.shop.poulten-graf.de



Glass Pasteur Pipettes



NEW:

Pasteur Pipettes
D810 + D812

made of
Borosilicate
Glass are now also
available.

Get information!



Superior
Quality!



volac Glass Pasteur Pipettes

All **volac** Glass Pasteur Pipettes are precision made from carefully selected tubing on fully automatic machines ensuring total uniformity of production. Tips and tapers are centred and can be checked by rolling the pipettes on a flat surface. Length and taper of the tips are carefully controlled. Constrictions for plugging are neatly formed at a uniform distance from the fire polished end. Constrictions are heat formed which avoids weaknesses or contamination of the glass found when moulding or tooling processes are used (ISO7712).

volac Glass Pasteur Pipettes are made with a slightly heavier wall thickness to provide the best quality product allied to greater strength to give lower breakages during transit and reduced risk of breakage in use or when fitting PVC teats. **volac** Glass Pasteur Pipettes are available either plain, pre-plugged with non-absorbent roving or pre-plugged and sterilised by vacuum autoclaving at 134 °C.

Plain and pre-plugged pasteurs are packed in inner boxes of 250 pieces with four boxes per standard outer pack of 1000 pieces.

Technical Details:

- available in standard sizes 150 mm, 230 mm, 270 mm, plain or plugged, sterilised or non-sterile
- made of soda lime glass or borosilicate glass (NEW)
- '50 dropper' pipettes also available
- superior wall thickness as standard
- pasteurs 230 mm long can be supplied with heavy wall
- all standard pipettes can be supplied with 'sealed points' and break marks on request

D810, D812, D 812 HW (length 230 mm)

	Standard Past. Pip. mm	Heavy wall Past.Pip. mm
Top Body outer ø	6.8 - 7.1	7.70 - 8.10
Wall thickness	0.5 - 0.53	0.80 - 0.90
Jet outer ø	1.50 ± 0.2	

volac Glass Disposable Pasteur Pipettes

	unplugged non-sterile soda lime glass	pre-plugged non-sterile soda lime glass
Length	Art.-No.	Art.-No.
150 mm	D810	D810/PP
230 mm	D812	D812/PP
270 mm	D814	D814/PP

	unplugged non-sterile borosilicate glass
Length	Art.-No.
150 mm	D810/BORO
230 mm	D812/BORO

NEW
made of
borosilicate glass

Packing unit: 1000 pcs
(4 inner packs x 250 pcs)

	pre-sterilised, pre-plugged soda lime glass
Length	Art.-No.
150 mm	D810/VA
230 mm	D812/VA

Packing unit: 250 pcs. (10 x 25)

